## **Dropping Low Phred Quality Score Sequences**

So in DNA Sequences, the Phred quality score is an important measure.

## **Phred Quality Score**

Base calling accuracy, measured by the Phred quality score (Q score), is the most common metric used to assess the accuracy of a sequencing platform.

It indicates the probability that a given base is called *incorrectly* by the sequencer.

## For example:

- If Phred assigns a Q score of 30 (Q30) to a base, this is equivalent to the probability of an incorrect base call 1 in 1000 times
  (Table 1). This means that the base call accuracy (i.e., the probability of a correct base call) is 99.9%.
- A lower base call accuracy of 99% (Q20) will have an incorrect base call probability of 1 in 100,
  meaning that every 100 bp sequencing read will likely contain an error.
- When sequencing quality reaches Q30, virtually all of the reads will be perfect, having zero errors and ambiguities. This is why Q30 is considered a benchmark for quality in next-generation sequencing.

Table 1: Quality Scores and Base Calling Accuracy

Phred Quality Score	Probability of Incorrect Base Call	Base Call Accuracy
10	1 in 10	90%
20	1 in 100	99%
30	1 in 1,000	99.9%
40	1 in 10,000	99.99%
50	1 in 100,000	99.999%

Reference: Illumina ~ Q-Scores

Since the phred quality score is for each base, even if a single base has a score < 30, we have to drop that part of the sequence. (part of a parquet batch file, 60.3 M reads were split into batches, so we got 603 batches of 1,00,000 sequences)

Results from the code:



Almost half the sequences got dropped

## **Next Steps**

- 1. Verify if this is the right approach
- 2. Loop this over all the batches to have the best quality sequences only (i.e 100 bp only and phred >= 30)
- 3. Go to the embeddings codes (DNABERT and DNABERT 2, done in the report prior to this)
- 4. Loop that over all the batches to have the embeddings
- 5. Compare embeddings quality for DNN training